


Product Datasheet

M6PR antibody (orb330308)

Description	Rabbit polyclonal antibody to M6PR
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human M6PR
Target	M6PR
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: RLKPLFNKSFESTVGQGSPTYIYIFRVCREAGNHTSGAGLVQINKSNGKE
Clonality	Polyclonal
MW	28kDa
Uniprot ID	P20645
NCBI	NP_002346
Expiration Date	12 months from date of receipt.



Western blot analysis of human, mouse ti...



Western blot analysis of Jurkat cell lys...