



Product Datasheet

ABCC8 antibody (orb330321)

Description	Rabbit polyclonal antibody to ABCC8
Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Porcine, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human ABCC8
Target	ABCC8
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: PLAFCGSENHSAAYRVDQGVLLNNGCFVDALNVVPHVFLLFITFPILFIGW
Clonality	Polyclonal
MW	177kDa
Uniprot ID	Q09428
NCBI	NP_000343
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T cell lysat...



Immunohistochemical staining of human Sm...