



Product Datasheet

Alkaline Phosphatase antibody (orb33033)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	ALPI
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	56812
Uniprot ID	P09923
NCBI	NP_001622.2
Dilution Range	WB: 1:1000, IHC-P: 1:100, IF/ICC: 1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human sm...



Immunohistochemical staining of paraffin...