



Product Datasheet

SLC39A5 antibody (orb330369)

Description	Rabbit polyclonal antibody to SLC39A5
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human SLC39A5
Target	SLC39A5
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: HYLAQLFGLYGENGTLTAGGLARLLHSLGLGRVQGLRLGQHGPLTGRAAS
Clonality	Polyclonal
MW	56kDa
Uniprot ID	Q6ZMH5
NCBI	NP_775867
Expiration Date	12 months from date of receipt.



Western blot analysis of human Placenta ...



Immunohistochemical staining of human Placenta ...