


Product Datasheet

TGFBI antibody (orb330378)

Description	Rabbit polyclonal antibody to TGFBI
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human TGFBI
Target	TGFBI
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: LKNNVSVNKEPVAEPDIMATNGVVHVITNVLPANRPQERGDELADSA
Clonality	Polyclonal
MW	46kDa
Uniprot ID	Q15582
NCBI	AAH26352
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat tissue u...



Western blot analysis of Jurkat cell lys...