



Product Datasheet

TCEB1 antibody (orb33038)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	ELOC
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	12473
Uniprot ID	Q15369
NCBI	NP_001191786.1 , NP_001191790.1 , NP_001191788.1 , NP_001191791.1 , NP_001191793.1 , NP_001191792.1 , NP_001191787.1 , NP_005639.1 , NP_001191789.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of TCEB1 antibody ...