

Product Datasheet

GPC3 antibody (orb33048)

Description	Rabbit polyclonal antibody to GPC3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	GPC3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	65563
Uniprot ID	P51654
NCBI	NP_004475.1 , NP_001158091.1 , NP_001158089.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

