



Product Datasheet

CPT1B antibody (orb330486)


Description	Rabbit polyclonal antibody to CPT1B
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat, Sheep, Zebrafish
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human CPT1B
Target	CPT1B
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: DLEMQFQRILDDPSPQPGEELKLAALTAGGRVEWAQARQAFFSSGKNKAA
Clonality	Polyclonal
MW	88 kDa
Uniprot ID	Q92523
NCBI	NP_004368
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human CPT1B



Western blot analysis of human CPT1B



Western blot analysis of HT1080 cell lysate