




Product Datasheet

CD9 antibody (orb33059)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	CD9
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	25416
Uniprot ID	P21926
NCBI	NP_001760.1
Dilution Range	IF/ICC: 1:25, WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of Human Cer...



Western blot analysis of in HepG2 cell l...