

Product Datasheet

SPINT2 antibody (orb330593)

Description	Rabbit polyclonal antibody to SPINT2
Species/Host	Rabbit
Reactivity	Human, Yeast
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human SPINT2
Target	SPINT2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: MLRCFRQQENPPLPLGSKVVVLGAGLFVMVLILFLGASMVYLIRVARRNQE
Clonality	Polyclonal
MW	28kDa
Uniprot ID	O43291
NCBI	NP_066925
Expiration Date	12 months from date of receipt.



Western blot analysis of human heart tis...