



Product Datasheet

SRPRB antibody (orb330594)


Description	Rabbit polyclonal antibody to SRPRB
Species/Host	Rabbit
Reactivity	Canine, Equine, Human, Mouse, Rat, Yeast, Zebrafish
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human SRPRB
Target	SRPRB
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: QDIAMAKSAKLIQQLEKELNLTNRVTRSAAPSTLYSSSTAPAQLGKKGKE
Clonality	Polyclonal
MW	30 kDa
Uniprot ID	Q9Y5M8
NCBI	NP_067026
Expiration Date	12 months from date of receipt.



Western blot analysis of human Jurkat ti...



Western blot analysis of HeLa cell lysat...



Immunohistochemical staining of human Br...