




Product Datasheet

LALBA antibody (orb33064)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	LALBA
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	16225
Uniprot ID	P00709
NCBI	NP_002280.1
Dilution Range	IHC-P: 1:10-50, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of CEM cell line lysate.



Western blot analysis of human breast tissue lysate.