

Product Datasheet

SAE1 antibody (orb330657)

Description	Rabbit polyclonal antibody to SAE1
Species/Host	Rabbit
Reactivity	Canine, Guinea pig, Human, Mouse, Rat, Yeast, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human SAE1
Target	SAE1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: MVEKEEAGGGISEEEAAQYDRQIRLWGLEAQKRLRASRVLLVGLKGLGAE
Clonality	Polyclonal
MW	38kDa
Uniprot ID	Q9UBE0
NCBI	NP_005491
Expiration Date	12 months from date of receipt.



Western blot analysis of human brain tis...