


Product Datasheet

ACTR1B antibody (orb330662)


Description	Rabbit polyclonal antibody to ACTR1B
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human ACTR1B
Target	ACTR1B
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: KIKISAPQERLYSTWIGGSILASLDTFKKMWVSKKEYEEDGSRAIHRKTF
Clonality	Polyclonal
MW	42kDa
Uniprot ID	P42025
NCBI	NP_005726
Expiration Date	12 months from date of receipt.

ACTR1B




Western blot analysis of human Fetal Liv...

ACTR1B



Western blot analysis of human Fetal Lun...

ACTR1B



Western blot analysis of human 721_B tis...