



Product Datasheet

TRIB1 antibody (orb330687)


Description	Rabbit polyclonal antibody to TRIB1
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat, Yeast
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human TRIB1
Target	TRIB1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: TEERTQLRLESLEDTHIMKGEDDALSDKHGCPAYVSPEILNTTGTYSGKA
Clonality	Polyclonal
MW	41 kDa
Uniprot ID	Q96RU8
NCBI	NP_079471
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Lu...



Western blot analysis of human COLO205 t...



Western blot analysis of human Intestine...