



Product Datasheet

ATG5 antibody (orb330734)


Description	Rabbit polyclonal antibody to ATG5
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human ATG5
Target	ATG5
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: DPEDGEKKNQVMIHGIEPMLETPLQWLSEHLSYPDNFLHISIIPQPTD
Clonality	Polyclonal
MW	32kDa
Uniprot ID	Q9H1Y0
NCBI	NP_004840
Expiration Date	12 months from date of receipt.



Western blot analysis of human Fetal Hepatoma (HeLa) cells.



Western blot analysis of HepG2 cell lysate.



Immunohistochemical staining of human placental tissue (PI).