


Product Datasheet

LYN antibody (orb330772)


Description	Rabbit polyclonal antibody to LYN
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human LYN
Target	LYN
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -20C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P07948
Expiration Date	12 months from date of receipt.



Western blot analysis of human Fetal Lysat...



Western blot analysis of Hela cell lysat...



Western blot analysis of human Fetal Bra...