


Product Datasheet


CYP26B1 antibody (orb330920)

Description	Rabbit polyclonal antibody to CYP26B1
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human CYP26B1
Target	CYP26B1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: LRWAATRDKSKLPIPKGSMGFPLIGETGHWLLQGSGFQSSRREKYGNVF
Clonality	Polyclonal
MW	57kDa
Uniprot ID	Q9NR63
NCBI	NP_063938
Expiration Date	12 months from date of receipt.

CYP26B1



Western blot analysis of human Fetal Bra...



Western blot analysis of Jurkat cell lys...