

Product Datasheet

ERCC4 antibody (orb330965)

Description	Rabbit polyclonal antibody to ERCC4
Species/Host	Rabbit
Reactivity	Canine, Equine, Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human ERCC4
Target	ERCC4
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: FLLRLYRKTFEKDSKAEVVMKFRKEDSSKRIRKSHKRPKDPQNKERAST
Clonality	Polyclonal
MW	101kDa
Uniprot ID	Q92889
NCBI	NP_005227
Expiration Date	12 months from date of receipt.



Western blot analysis of PANC1 cell lysa...