



Product Datasheet

ALDH9A1 antibody (orb331054)


Description	Rabbit polyclonal antibody to ALDH9A1
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human ALDH9A1
Target	ALDH9A1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: MGPLINRPHLERVLGFVKVAKQGAQVLCGGDIYVPEDPKLKDGYYMRPC
Clonality	Polyclonal
MW	56kDa
Uniprot ID	B9EKV4
NCBI	NP_000687
Expiration Date	12 months from date of receipt.



Western blot analysis of human Liver tis...



Western blot analysis of human Fetal Bra...



Western blot analysis of human Placenta ...