

Product Datasheet

KCTD18 antibody (orb33113)

Description	KCTD18 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human BTB/POZ domain-containing protein KCTD18 protein (1-200AA)
Target	KCTD18
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6PI47
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.