


Product Datasheet

ORC3 antibody (orb331176)

Description	Rabbit polyclonal antibody to ORC3 also known as Origin recognition complex subunit Latheo, whose ex...
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Porcine, Rat, Yeast, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Target	ORC3
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 - 1 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: VVTAAEKMDANSATSEEMNEIIHARFIRAVSELELLGFIKPTKQKTDHVA
Clonality	Polyclonal
MW	82kDa
Uniprot ID	Q9UBD5
NCBI	NP_862820
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 Whole Cell Lysate