

Product Datasheet

SFTPA1 antibody (orb331250)

Description	Rabbit polyclonal antibody to SFTPA1
Species/Host	Rabbit
Reactivity	Canine, Equine, Human, Mouse, Porcine, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of human SFTPA2
Target	SFTPA2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: KYNTYAYVGLTEGSPGDFRYSDGTPVNYTNWYRGEPAGRGKEQCVEMYT
Clonality	Polyclonal
MW	20kDa
Uniprot ID	Q8IWL1
NCBI	NP_001092138.1
Expiration Date	12 months from date of receipt.



Western blot analysis of human HCT15 cells...