

Product Datasheet

LYPLAL1 antibody (orb331336)

Description	Rabbit polyclonal antibody to LYPLAL1
Species/Host	Rabbit
Reactivity	Canine, Equine, Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of Human LYPLAL1
Target	LYPLAL1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: MCQVLTDLIDEEVKSGIKKNRILIGGFSMGGCMAMHLAYRNHQDVAGVFA
Clonality	Polyclonal
MW	26kDa
Uniprot ID	Q5VWZ2
NCBI	NP_620149
Expiration Date	12 months from date of receipt.



Western blot analysis of COLO205 Whole Cell Lysate