

Product Datasheet

SNX15 antibody (orb331339)

Description	Rabbit polyclonal antibody to SNX15
Species/Host	Rabbit
Reactivity	Canine, Equine, Goat, Guinea pig, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N-terminal region of Human SNX15
Target	SNX15
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: RKLHGDLAYTHRNLFRRLEEFPAFPRAQVGRFEASVIEERRKGAEDLLR
Clonality	Polyclonal
MW	29kDa
Uniprot ID	E5KQS6
NCBI	NP_680086
Expiration Date	12 months from date of receipt.



Western blot analysis of THP-1 Whole Cell...