


Product Datasheet

AP1B1 antibody (orb331499)

Description	Rabbit polyclonal antibody to AP1B1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of Human AP1B1
Target	AP1B1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: PSAFVEGGRGVHKSLPPRTASSESAESPETAPTGAPPGEQPDVIPAQGD
Clonality	Polyclonal
MW	104kDa
Uniprot ID	Q10567
Expiration Date	12 months from date of receipt.



Western blot analysis of human 293T cells. The blot shows a single band at approximately 104 kDa, corresponding to the AP1B1 protein. Molecular weight markers are indicated on the left at 168 kDa, 144 kDa, 90 kDa, 65 kDa, and 40 kDa.