

Product Datasheet

PKP3 antibody (orb33152)

Description	Rabbit polyclonal antibody to PKP3. The pretein is a member of the arm-repeat & plakophilin gene fam...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PKP3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	87082
Uniprot ID	Q9Y446
NCBI	NP_009114.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

