



Product Datasheet

FXVD5 antibody (orb33154)

Description	Rabbit polyclonal antibody to FXVD5. The protein belongs to the FXVD family, is involved in down-reg...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	FXVD5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	19472
Uniprot ID	Q96DB9
NCBI	NP_001158077.1 , NP_659003.1 , NP_054883.3
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of CEM cell line l...