



Product Datasheet

Ptk6 antibody (orb33171)

Description	Rabbit polyclonal antibody to Ptk6. The protein, is also known as Breast tumor Kinase (BRK), is a Sr...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	Ptk6
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	51972
Uniprot ID	Q64434
NCBI	NP_033210.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse heart tis...