

Product Datasheet

Mouse Hoxc12 antibody (orb33173)

Description Rabbit polyclonal antibody to Hoxc12. The protein belongs to the Abd-B homeobox family. The homeobox...

Species/Host Rabbit

Reactivity Mouse

Conjugation Unconjugated

Tested Applications WB

Immunogen Synthetic Peptide

Target Hoxc12

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 30345

Uniprot ID [Q8K5B8](#)

NCBI [NP_034593.1](#)

Dilution Range WB: 1:1000

Expiration Date 12 months from date of receipt.

