




Product Datasheet

ALPP antibody (orb33185)

Description	Rabbit polyclonal antibody to ALPP. The protein belongs to the alkaline phosphatase family, and is l...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	ALPP
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	57954
Uniprot ID	P05187
NCBI	NP_001623.3
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cell line...



Western blot analysis of 293 cell lysate...