




Product Datasheet

SIGLEC6 antibody (orb33275)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	SIGLEC6
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	49913
Uniprot ID	O43699
NCBI	NP_001171018.1 , NP_001171019.1 , NP_942143.3 , NP_001171020.1 , NP_942142.3 , NP_001236.4
Dilution Range	IHC-P: 1:50-100, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 cell line l...



Western blot analysis of in mouse spleen...