



Product Datasheet

Estrogen Receptor beta antibody (orb33298)

Description	Rabbit polyclonal antibody to Estrogen Receptor beta
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	ESR2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	59216
Uniprot ID	Q92731
NCBI	NP_001258805.1 , NP_001035365.1 , NP_001258806.1 , NP_001428.1 , NP_001201832.1 , NP_001201831.1
Dilution Range	WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
HepG2 cells
...