



Product Datasheet

LC3 antibody (orb33318)

Description	Rabbit polyclonal antibody to LC3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant Protein
Target	This LC3 antibody is generated from rabbits immunized with full-length recombinant human LC3 (APG8a).
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q9H492 , Q9GZQ8
NCBI	NP_115903.1 , NP_852610.1
Dilution Range	IF/ICC: 1:200, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 cell lysate...