


Product Datasheet

NAPEPLD antibody (orb33326)

Description	NAPEPLD antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human N-acyl-phosphatidylethanolamine-hydrolyzing phospholipase D protein (1-200AA)
Target	NAPEPLD
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6IQ20
Dilution Range	WB:1:1000-5000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
Mouse large
in...