




Product Datasheet

Cleaved LC3 antibody (orb33336)

Description	Rabbit polyclonal antibody to Cleaved LC3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	This Cleaved LC3A antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 89-120 amino acids from human Cleaved LC3A or LC3B.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q9H492, Q9GZQ8
NCBI	NP_115903.1, NP_852610.1
Dilution Range	IF/ICC: 1:25, IHC-P: 1:10-50, WB: 1:1000, IF/ICC-1:10-50
Expiration Date	12 months from date of receipt.



WB analysis of A431 cell line using Clea...



WB analysis of A431 cell line using Clea...