




Product Datasheet

Cleaved LC3 antibody (orb33337)

Description	Rabbit polyclonal antibody to Cleaved LC3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Synthetic Peptide
Target	MAP1LC3B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	14688
Uniprot ID	Q9GZQ8
NCBI	NP_073729.1
Dilution Range	IF/ICC: 1:10-50, IF/ICC-1:10-50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of SY5Y cell...



Immunofluorescent analysis of mouse leukocytes...