

Product Datasheet

Cytokeratin 12 antibody (orb333741)

Description	Rabbit polyclonal antibody to Cytokeratin 12
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N-terminal region of Human KRT12
Target	KRT12
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -20C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q99456
Expiration Date	12 months from date of receipt.



Western blot analysis of human MCF7 Whole Cell Lysate (WCL) probed with Cytokeratin 12 antibody (orb333741).