



Product Datasheet

Tlr2 Antibody (orb333869)


Description	Mouse monoclonal antibody to Tlr2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC
Immunogen	Purified recombinant fragment of mouse Tlr2 (AA: 628-777) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	90kDa
Entrez	24088
Dilution Range	ICC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse recombinant...



Western blot analysis of HEK293 using Tl...



ELISA Black line: Control Antigen (100 n...