



Product Datasheet

RNLS Antibody (orb333878)


Description	Mouse monoclonal antibody to RNLS
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified recombinant fragment of human RNLS (AA: 68-242) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	37.8kDa
Entrez	55328
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 using RN...



ELISA Black line: Control Antigen (100 n...



Western blot analysis of human recombina...