


Product Datasheet


RF1 Antibody (orb333891)

Description	Mouse monoclonal antibody to RF1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Immunogen	Purified recombinant fragment of human RF1 (AA: 288-437) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	49kDa
Entrez	2107
Dilution Range	WB: 1/500 - 1/2000; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.




Flow
cytometric
analysis of
HepG2 cells

...



Western blot
analysis of
human
recombina...



Western blot
analysis of
HEK293 using
RF...