



Product Datasheet

TWIST2 Antibody (orb333899)


Description	Mouse monoclonal antibody to TWIST2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified recombinant fragment of human TWIST2 (AA: 1-160) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	18.1kDa
Entrez	117581
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 using TWIST2 antibody.



ELISA Black line: Control Antigen (100 ng) n...



Western blot analysis of human recombinant TWIST2 protein.