

Product Datasheet

NLRC4 (phospho-Ser-533) (orb333906)

Description Mouse monoclonal antibody to NLRC4 (phospho-Ser-533)

Species/Host Mouse

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications ELISA, ICC, WB

Immunogen Synthesized peptide of mouse phospho-NLRC4(Ser-533) (AA: 525-538) expressed in E. Coli.

Preservatives Purified antibody in PBS with 0.05% sodium azide

Concentration 1mg/ml

Storage 4°C; -20°C for long term storage. Avoid freeze/thaw cycles.

Note For research use only

Isotype Mouse IgG2b


Clonality Monoclonal

MW 116.7kDa


Entrez [268973](#)

Dilution Range WB: 1/500 - 1/2000; ICC: 1/200 - 1/1000; ELISA: 1/10000


Expiration Date 12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Western blot analysis of NIH/3T3 (Lane1)...



ELISA Black line: Control Antigen (100 ng)