



Product Datasheet

PRDM4 Antibody (orb333920)


Description	Mouse monoclonal antibody to PRDM4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Purified recombinant fragment of human PRDM4 (AA: 476-575) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	88kDa
Entrez	11108
Dilution Range	RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of human recombinant PRDM4.



Western blot analysis of HEK293 cells using PRDM4 antibody.



ELISA Black line: Control Antigen (100 ng).