



## Product Datasheet

### LRP1 Antibody (orb333966)


|                            |                                                                                |
|----------------------------|--------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to LRP1                                              |
| <b>Species/Host</b>        | Mouse                                                                          |
| <b>Reactivity</b>          | Human                                                                          |
| <b>Conjugation</b>         | Unconjugated                                                                   |
| <b>Tested Applications</b> | ELISA, ICC, IHC                                                                |
| <b>Immunogen</b>           | Purified recombinant fragment of human LRP1 (AA: 20-155) expressed in E. Coli. |
| <b>Preservatives</b>       | Purified antibody in PBS with 0.05% sodium azide                               |
| <b>Concentration</b>       | 1mg/ml                                                                         |
| <b>Storage</b>             | 4°C; -20°C for long term storage. Avoid freeze/thaw cycles.                    |
| <b>Note</b>                | For research use only                                                          |
| <b>Isotype</b>             | Mouse IgG1                                                                     |
| <b>Clonality</b>           | Monoclonal                                                                     |
| <b>MW</b>                  | 504kDa                                                                         |
| <b>Entrez</b>              | <a href="#">4035</a>                                                           |
| <b>Dilution Range</b>      | IHC: 1/200 - 1/1000; ICC: 1/200 - 1/1000; ELISA: 1/10000                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                |



Immunohistochemical analysis of formalin...



Immunohistochemical analysis of formalin...



Western blot analysis of human recombina...