



Product Datasheet

PRDM1 Antibody (orb333970)


Description	Mouse monoclonal antibody to PRDM1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC
Immunogen	Purified recombinant fragment of human PRDM1 (AA: 690-825) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG2a
Clonality	Monoclonal
MW	91.7kDa
Entrez	639
Dilution Range	ICC: 1/100 - 1/500; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of human recombinant PRDM1.



Western blot analysis of HEK293 cells using PRDM1 antibody.



ELISA Black line: Control Antigen (100 ng/ml).