



Product Datasheet

PDGFRB Antibody (orb333998)


Description	Mouse monoclonal antibody to PDGFRB
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Purified recombinant fragment of human PDGFRB (AA: 33-133) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	124kDa
Entrez	5159
Dilution Range	RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of human recombinant...



Western blot analysis of HEK293 using PD...



ELISA Black line: Control Antigen (100 n...