



Product Datasheet

KDM5B Antibody (orb334007)


Description	Mouse monoclonal antibody to KDM5B
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Purified recombinant fragment of human KDM5B (AA: 231-319) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	175.6kDa
Entrez	10765
Dilution Range	RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of human recombinant KDM5B.



Western blot analysis of HEK293 cells using KDM5B antibody.



ELISA Black line: Control Antigen (100 nM).