



Product Datasheet

Neuropilin-1 Antibody (orb334026)


Description	Mouse monoclonal antibody to Neuropilin-1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Immunogen	Synthesized peptide of human Neuropilin-1 (AA: 45-59).
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	103.1kDa
Entrez	8829
Dilution Range	WB: 1/500 - 1/2000; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat (Lane1),...



ELISA Black line: Control
Antigen (100 n...



Flow cytometric analysis of Hela cells u...