



Product Datasheet

THAP1 Antibody (orb334030)


Description	Mouse monoclonal antibody to THAP1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Purified recombinant fragment of human THAP1 (AA: 1-213) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	24.9kDa
Entrez	55145
Dilution Range	RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of human recombinant THAP1.



Western blot analysis of HEK293 cells using THAP1 antibody.



ELISA Black line: Control Antigen (100 ng).